




Vision Monitor

MonBaby^{Box}

The **MonBaby^{Box}** is a stand-alone system for recording flash ERGs and VEPs.



Mechanical specifications

1/1

For one optical head:

Width	400 mm
Height	200 mm
Depth	370 mm
Weight	9.5 kg

Electrical specifications

230 V, 0.7A or 110V, 1.4A 50 or 60 Hz

The **MonBaby^{Box}** contains its own isolation transformer.

It includes one power-outlet that can be used to power a laptop.

It can also be connected to a standard PC tower with a color monitor but, in that case, it requires an external isolation transformer.



WARNING: Do not use it in the presence of flammable anesthetics.